

# Eliminate

Complete disinfection in 4 minutes!



## About this product

Eliminate is a one-step disinfectant wipe, with no need to pre-wash or rinse the surface. It handles daily cleaning and disinfection without requiring additional sprays, sponges or buckets. Eliminates 99.9% of contaminants in just 15 seconds and complete disinfection in only 4 minutes; giving you the peace of mind that your surfaces are free of bacteria, viruses and fungi. The fast and efficient formula makes it ideal for busy environment where there is no time to waste on slow disinfecting time. The textured fabric of the wipe makes it ideal to absorb the chemistry, while capturing the dirt, bacteria and viruses without smearing them. Easy to use, its safe formula can be used on a wide variety of non-porous surfaces in both health care and non-health care settings.

**DIN: 02459507**

- Delivers effective cleaning, sanitizing, disinfecting and malodour control
- Ideal for use by non-custodial workers
- Worry free disinfecting
- Convenient and easy to use, saves considerable time cleaning
- Eliminates the risk of dilution errors – every single wipe contains the proper
- Active ingredients to ensure complete disinfection
- Eliminates the risk of cross-contamination from the conventional bucket method
- Disposable to prevent the spread of microorganisms and infections
- No harsh bleach scent
- Can be used as a floor cleaning/disinfectant wipe



### Kill Claims

Eliminate Wipes are a one-step disinfectant and bactericidal according to the current AOAC Use-Dilution Test Method modified in the presence of **5% organic serum against:**

#### Bacteria:

Staphylococcus aureus, Salmonella enterica, Pseudomonas aeruginosa

Escherichia coli O157:H7, ESBL Escherichia coli, Klebsiella pneumoniae, Listeria monocytogenes, Streptococcus pyogenes

#### Viruses:

HIV-1, Influenza Type A Virus / Brazil Flu, Influenza A Virus H1N1, 2013 Influenza A Virus (H7N9), Human HBV, Human HCV

#### Fungi:

Candida albicans



Code:

53377 (180 wipes)  
\*6 containers per box

